## **Lumin Test**

### **Microorganism contamination measurement kit**

*Catalogue No. BLT-0025*  
*Packing specification: 25 kits*

#### **Intended Use**

Lumin Test is an ATP test kit for microorganism assessment of surface. Lumin Test is a user-friendly kit use for cleaning verification of surfaces on processing equipment and other environments. Lumin Test can be used in a wide range of industries such as, food processing industry or services, health care, hospitality, sterile services in medical, and so on.

#### **Test Principle**

Adenosine Triphosphate (ATP) is the energy molecule found in all living cells. Measuring the amount of ATP indicates the amount of living cells on a surface and therefore can determine if a surface is clean or not. A special enzymatic reaction is introduced into the ATP, and then causes light emission. This emitted light is measured by a portable luminometer and is proportional to the amount of living organisms within the sample. The luminometer reports results in Relative Light Units (RLU). RLU result provides information on the level of infection within 15 seconds. The higher the RLU value, the more ATP present, and the dirtier the surface.

#### **Storage & Shelf Life**

- Store at 2 – 8 °C (36 – 46 °F) and out of direct sunlight
- Lumin Test will tolerate room temperature (< 25 °C) for 4 weeks
- Lumin Test have a 15-month shelf life

#### **Safety & Precautions**

- Components of Lumin Test do not pose any health risk
- Do not reuse a Lumin Test

**Kit components:** 25 swab test in aluminum pouch and instruction sheet.


#### **Test procedure**

Before beginning, turn on luminometer and Allow Lumin Test to equilibrate to room temperature (20- 25°C).

1- Pull top of swab out of swab tube.

2. Thoroughly swab a standard 10 x 10 cm (4 x 4 inches) area. Note: Do not touch swab or inside of sample device with fingers. Rotate swab vertically, horizontally while collecting sample to maximize sample collection on swab tip. Apply sufficient pressure to create flex in swab shaft.

3. After swabbing, replace swab back in swab tube.
4. Hold swab tube firmly and use thumb and forefinger to break upper part of red colored pieces by bending accordion-like part forward and backward. Crush accordion-like part twice, expelling all liquid down swab shaft.
5. Rinse swab tube in liquid by shaking for 5 –10 seconds.
- 6- Read sample in luminometer within 15 seconds. Note: Holding luminometer upright, insert entire Luminn Test into luminometer.

### Interpretation of Results

Scale of clearness level based on RLU measured is listed at the below table.

ATP Concentration (RLU)	Clearness level
< 10	Good
11-30	Average
> 30	Bad

#### Technical Support

Tel: 021-66428602; 66944047-49

Fax: 021- 66593861

Email: baharzist@yahoo.com

No. 9, 3th Golshid, Nakhlestan Blvd, Baharestan Blvd, Shams Abad Industrial Zone, Tehran, Iran

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