

TLC Kit


Multiple Drugs of Abuse Toxicology Kit

Cat No. BTK-52131

TLC kit is designed to simultaneously detect **8** different drugs via thin layer chromatography method in urine specimen of persons. TLC kit is used for screening and confirmation of drugs of abuse.

Urine Cut-Off Levels

Cut Off	Drugs
300 ng/ml	Morphine, Codeine, Methadone, Buprenorphine Tramadol
500 ng/ml	Amphetamine, Methamphetamine, Ecstasy



Chromatographic Patters

Advantages

- Fast and easy to use
- One TLC test for detection of **8** drugs of abuse, including Morphine, Codeine, Methadone, Buprenorphine, Tramadol, Amphetamine, Methamphetamine, Ecstasy
- Respose in **50** minutes.
- Qualitative analysis by chromatographic patterns of spot for each drugs
- Quantitative analysis with cut off values are as above table

