

Lumin Test

Cat No. BLT-0025

Microorganism contamination measurement kit

Lumin Test is an ATP (adenosine triphosphate) test kit for microorganism assessment of a surface. Lumin Test is used for hygiene monitoring, providing results in **15** seconds of inserting the sample into the luminometer device. The use of small pen sized sample swab makes taking the sample as easier and time efficient as ever.

ATP Concentration (RLU)	Clearness level
< 10	Good
11-30	Average
> 30	Bad



Applications

- Clean verification
- Measure the effectiveness of disinfection
- Immediate feedback on equipment or operators cleanliness

Customers

- Health care, hospitality and sterile services in medical
- Food process and production industry such as, dairy, fruit juices, canned food, meat, poultry, and seafood
- Food services such as, bakery, restaurants, cafeterias, sandwich shops coffee shops, fast food retail

